

Wireless *MOSkin* Dosimetry for In-Vivo Rectal Monitoring During HDR Prostate Brachytherapy

D. Cutajar^{1,2}, I. Peters¹, R. Brown², A. Howie², J. Poder^{1,2},
J. Bucci², A. Pogossov¹, A. Rosenfeld¹

1. Centre for Medical Radiation Physics, School of Physics, University of Wollongong, Australia
2. St. George Cancer Care Centre, St. George Hospital, Kogarah, Australia

CENTRE FOR
MEDICAL RADIATION PHYSICS



UNIVERSITY OF
WOLLONGONG






The MOSkin Dosimeter



- MOSkin™ detector is a novel MOSFET detector prototyped and developed at the Center for Medical Radiation Physics (CMRP), University of Wollongong, Australia.

- **The MOSkin™ advantages:**



- Small physical size, SV less than 1um
- Can be placed against cavity interfaces
- Real time read out
- Can be calibrated using Strontium-90
- Disposable after use



MOSkin™
Detector

Previous In-Vivo Rectal Measurements (St. G)

- 2007, 2008 – Performed using a prototype MOSkin within a sterilisable rubber tube, glued to a MEDRAD rectal balloon (Telisis 5 medical grade adhesive), placed within a condom



Previous In-Vivo Rectal Measurements (St. G)




- 2010 – Performed using *MOSkin* dosimeters mounted within RadiaDyne rectal balloons during manufacture



- Both trials showed promise, however the accuracy of measurement was limited by the movement of balloons in the rectum
 - Shifting, deflation

Dual Purpose Probe (DPP)

- Developed in collaboration with CMRP, SGCCC and the IRCSS (Milan, Italy)




M. Cararra, et al, Radiotherapy and Oncology, 118 (2016) 148-153

- MOSkins mounted on TRUS probe (rigid)
- Constant reference frame between imager (TRUS) and detectors
- Used for HDR prostate brachytherapy and gynaecological procedures (thin rigid probe)

DPP Clinical Use




- In Milan, Italy
- Real time TRUS based HDR prostate brachytherapy
- 18 fractions (36 measurement points)
- The measured doses were compared to
 - Pre-treatment image plan ($-2.1\% \pm 8.3\%$)
 - Post-treatment image plan ($-0.6 \pm 4.1\%$)
- This study showed the delivery needs to be performed within as short a time as possible



Similarly for Gynaecological

- Rectal wall doses were measured for a total of 22 measurements in 9 vaginal HDR brachytherapy fractions, for a total of 3 patients
- 3 detectors were attached to rectal catheter



MOSkin pigtails,
plugged into
reader

Orientation marker

Probe filled with gel

Radiopaque marker

3 MOSkin dosimeters
Kapton tails allow for
dosimeters to be wrapped
around catheter

New MOSkin Reader


- Improved measurement accuracy
 - 0.125 mV minimum signal as opposed to 1 mV
- Short cable connection from MOSkin to reader
 - Up to 6 devices at once
 - Battery powered (rechargeable)
- Wireless connection from reader to computer
 - Custom reader software
 - Bluetooth or 434 MHz connection
- Automatic temperature fluctuation compensation
 - Constantly measures leakage current through an on-chip diode
 - Offsets all measurements to account for temperature fluctuations



Temperature Calibration




- Performed for every *MOSkin* probe
- Probe placed on an inverted Peltier cooler (to heat it)
- As device cools, a trace of leakage current vs threshold voltage (*MOSkin* signal) is obtained
- The gradient of this is the temperature calibration
- By constantly comparing the leakage current to the threshold voltage, any changes in threshold voltage due to temperature fluctuation may be separated from the signal induced by radiation



Phantom Trial

- A DPP was created with three MOSkin probes attached to a TRUS
- Designed to span 5cm (base, mid-gland, apex)
- A CIRS prostate phantom was implanted under TRUS guidance
- Real-time TRUS based planning used
- Treatment delivered with MOSkin measurements obtained with the classic reader and the wireless reader



Phantom Trial



Position	Classic Reader	Wireless Reader
Base (0 mm)	0.6%	3.7%
Mid-Gland (24 mm)	4.3%	3.0%
Apex (50 mm)	2.6%	3.4%

Discrepancies from treatment plan


It was established that the wireless reader would supplant the classic reader with minimal differences in accuracy




Clinical Trials



- Real-time TRUS based planning for HDR Prostate BT introduced at StG CCC Late 2016
 - 2 Fractions, 1 week apart
- MOSkins to be used to verify rectal doses
- 4 MOSkins mounted to the side of the TRUS probe (1.5cm pitch)
- Probe placed inside TRUS balloon (imaging, hygiene)
- Trials to begin late Feb 2018
- 13 patients, 26 fractions, 104 detectors (6 months)



Conclusions



- Development of the MOSkin system has progressed to improve
 - Accuracy (better electronics)
 - Useability (wireless, software based acquisition)
 - Stability (temperature correction, lower operating voltage)
- Clinical trials will proceed to
 - Incorporate 4 MOSkins onto a TRUS
 - Measure the rectal wall dose during TRUS based HDR BT
 - Physically validate the real-time TRUS based HDR BT procedure at St.G CCC

